



Product Datasheet

PPP2R4 Antibody (orb97413)


Description	Mouse monoclonal antibody to PPP2R4
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IHC
Immunogen	Purified recombinant fragment of human PPP2R4 expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	41kDa
Entrez	5524
Dilution Range	IHC: 1/200 - 1/1000; ICC: 1/200 - 1/1000; RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



IHC-P of lung cancer tissues using PPP2R...



Confocal immunofluorescence analysis of ...



Immunohistochemical analysis of formalin...