



Product Datasheet

PIK3CA (orb97439)


Description	Mouse monoclonal antibody to PIK3CA
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Entrez	5290
Expiration Date	12 months from date of receipt.



Confocal immunofluorescence analysis of ...



Western blot analysis of HEK293 (Lane1) ...



Flow cytometric analysis of Jurkat cells...