



Product Datasheet

RAP1A Antibody (orb97460)


Description	Mouse monoclonal antibody to RAP1A
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	Purified recombinant fragment of human RAP1A expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	43.2kDa
Entrez	5906
Dilution Range	RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of HeLa cells U...



Western blot analysis of HEK293 (Lane1) ...



Western blot analysis of human (AA: 28-1...