



Product Datasheet

INCENP Antibody (orb97493)


Description	Mouse monoclonal antibody to INCENP
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IHC
Immunogen	Purified recombinant fragment of human INCENP (AA: 369-583) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	105.4kDa
Entrez	3619
Dilution Range	IHC: 1/200 - 1/1000; ICC: 1/200 - 1/1000; RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



IHC-P of rectum cancer tissues using INCENP antibody



Confocal immunofluorescence analysis of INCENP antibody



Line graph illustrates about the Ag-Ab reaction