



Product Datasheet

MYL3 Antibody (orb97569)


Description	Mouse monoclonal antibody to MYL3
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Purified recombinant fragment of MYL3 expressed in E. Coli.
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	22kDa
Entrez	4634
Dilution Range	WB: 1/500 - 1/2000; IHC: 1/200 - 1/1000; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis of rat fetal heart...



Immunohistochemical analysis of formalin...



IHC-P of human Heart tissues using MYL3 ...