



Product Datasheet

IGFR(Phospho-Tyr1161) Antibody (orb97626)


| | |
|----------------------------|---|
| Description | IGFR(Phospho-Tyr1161) Antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | IF, IHC, WB |
| Storage | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| Note | For research use only |
| Clonality | Recombinant |
| Uniprot ID | P08069 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of IGF-1R (Ab-1161...



Immunohistochemical analysis of formalin...



Confocal immunofluorescence analysis of ...