



Product Datasheet

FKHRL1(Phospho-Ser253) Antibody (orb97629)


Description	FKHRL1(Phospho-Ser253) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	043524
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from N...



Immunohistochemical analysis of formalin...



Confocal immunofluorescence analysis of ...