




Product Datasheet

ALDH2 antibody (orb97662)

Description	Mouse monoclonal to ALDH2. This ALDH2 monoclonal antibody is generated from mouse immunized with ALDH...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Recombinant Protein
Target	ALDH2
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG1, Igk
Clonality	Monoclonal
MW	56381
Uniprot ID	P05091
NCBI	NP_001191818.1 , NP_000681.2
Dilution Range	WB: 1:1000, IHC-P: 1:50-100, FACS: 1:10-50, IF/ICC: 1:10-50
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HepG2 cel...



Immunohistochemical staining of paraffin...