




Product Datasheet

ALDH1A1 antibody (orb97664)

Description	Mouse monoclonal antibody to ALDH1A1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant Protein
Target	ALDH1A1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, eluted with high and low pH buffers and neutralized immediately, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	54862
Uniprot ID	P00352
NCBI	NP_000680.2
Dilution Range	WB: 1:500-1000, IF/ICC: 1:10-50, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HepG2 cel...



Immunofluorescence analysis of NCI-H460 ...