




Product Datasheet

CD44 antibody (orb97676)

Description	Mouse monoclonal to CD44.This CD44 monoclonal antibody is generated from mouse immunized with CD44 r...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Target	CD44
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG2a ,K
Clonality	Monoclonal
MW	81538
Uniprot ID	P16070
NCBI	NP_001001391.1 , NP_001189484.1 , NP_001189485.1 , NP_001001392.1 , NP_001001389.1 , NP_001189486.1 , NP_001001390.1 , NP_000601.3
Dilution Range	IHC-P: 1:10-50, WB: 1:2000, IF/ICC: 1:10-50, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Immunohistochemical staining of paraffin...