



Product Datasheet

SOX2 antibody (orb97681)


Description	Mouse monoclonal to SOX2. SOX2 recombinant protein is used to produce this monoclonal antibody.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Recombinant Protein
Target	SOX2
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	34310
Uniprot ID	P48431
NCBI	NP_003097.1
Dilution Range	WB: 1:200-2000, IF/ICC: 1:100, IHC-P: 1:50-100, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of SOX2 protein...



Western blot analysis of recombinant pro...



Flow cytometric analysis of NCI-H460 cel...