
Product Datasheet

Myoglobin/Cy3 antibody (orb976863)

Description

Mouse monoclonal antibody to

Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Cy3
Tested Applications	ICC, IF
Immunogen	Purified Human Myoglobin (Natural Human Myoglobin protein)
Target	MB
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	17kDa
Uniprot ID	P02144
Dilution Range	IF=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.