


## Product Datasheet

### hMSRB3 antibody (orb97722)

<b>Description</b>	Anti-human-MSRB3 (h-Methionine sulfoxide reductase B3) is a rabbit polyclonal antibody that was prod...
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	Unconjugated
<b>Form/Appearance</b>	liquid (Supplied as serum, preserved in glycerol)
<b>Storage</b>	store at -20 °C
<b>Note</b>	For research use only
<b>Application notes</b>	Anti-human MsrA can be used for Western blot, ELISA and Immunocytochemistry.
<b>Clonality</b>	Recombinant
<b>Expiration Date</b>	12 months from date of receipt.



Immunofluorescence analysis of using hMS...