

Product Datasheet

SHC (phospho-Tyr349) antibody (FITC) (orb9784)

Description	Rabbit polyclonal to SHC (phospho-Tyr349) antibody conjugated to FITC
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human SHC around the phosphorylation site of Tyr349 HQ(p-Y)YN
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	63kDa
Uniprot ID	P29353
Dilution Range	IF: 1:50-200
Expiration Date	12 months from date of receipt.