

Product Datasheet

SHP2 (phospho-Tyr580) antibody (FITC) (orb9796)

Description	Rabbit polyclonal antibody to Anti-Phospho-SHP2 (Tyr580)/FITC
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	FITC
Tested Applications	ICC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human PTPN11/SHP2 around the phosphorylation site of Tyr580 RV(p-Y)E
Target	PTPN11
Preservatives	Lyophilized or Liquid
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	66kDa
Uniprot ID	Q06124
Dilution Range	ICC/IF=1:25
Expiration Date	12 months from date of receipt.