

## Product Datasheet

# Alcohol dehydrogenase antibody (Biotin) (orb98141)

<b>Description</b>	Rabbit polyclonal to Alcohol dehydrogenase (Biotin)
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	Biotin
<b>Storage</b>	Storage: Lyophilized – can be stored at +4°C. After reconstitution aliquots can be stored at -2°C
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">P12762</a>
<b>Hazard Information</b>	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals. This datasheet is as accurate as reasonably achievable, but Nordic-MUBio accepts no liability for any inaccuracies or omissions in this information.
<b>Expiration Date</b>	12 months from date of receipt.