

Product Datasheet

Donkey anti Goat IgG (TRITC) antibody (orb98156)

Description	Tetramethylrhodamine (TMR) and its isothiocyanate derivative-conjugated IgG fraction of polyclonal d...
Species/Host	Donkey
Conjugation	TRITC
Storage	Storage: Lyophilized – can be stored at +4°C. After reconstitution aliquots can be stored at -2°C
Note	For research use only
Clonality	Polyclonal
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals. This datasheet is as accurate as reasonably achievable, but Nordic-MUbio accepts no liability for any inaccuracies or omissions in this information.
Expiration Date	12 months from date of receipt.