
Product Datasheet

Goat anti Human IgG IgA IgM antibody (orb98159)

Evolve Biosystems.
Description

Goat polyclonal to IgG IgA IgM

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

Immunofixation

Applications

Preservatives

Goat Ig lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity.

Form/Appearance

Goat Ig lyophilized from a solution in phosphate buffered saline (PBS, pH 7.2). No preservative added, as it may interfere with the antibody activity.

Storage

Storage: The lyophilized product is shipped at ambient temperature and may be stored at +4°C; prolonged storage at or below -20°C. Prior to use, an aliquot is thawed slowly at ambient temperature, spun down again and used to prepare working dilutions by adding sterile phosphate buffered saline (PBS, pH 7.2). Repeated thawing and freezing should be avoided. Working dilutions should be stored at +4°C, not refrozen, and preferably used the same day. If a slight precipitation occurs upon storage, this should be removed by centrifugation. It will not affect the performance of the product

Note

For research use only

Application notes

ELISA, Immunocytochemistry, Immunohistochemistry, (In)direct immunofluorescence.

Clonality

Polyclonal

Source

Highly purified IgG, and pools of homogenous IgA, IgM and IgD isolated from human serum. Freund's complete adjuvant is used in the first step of the immunization procedure.

Hazard Information

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals. This datasheet is as accurate as reasonably achievable, but Nordic-MUbio accepts no liability for any inaccuracies or omissions in this information.

Expiration Date

12 months from date of receipt.