

## Product Datasheet

### Goat Anti Rabbit IgG (orb98806)

|                          |   |
|--------------------------|---|
| <b>Description</b>       | Goat Anti Rabbit IgG  |
| <b>Species/Host</b>      | Goat  |
| <b>Conjugation</b>       | Unconjugated  |
| <b>Form/Appearance</b>   | Purified  |
| <b>Storage</b>           | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>              | For research use only   |
| <b>Application notes</b> | Suggested applications: EIA, CLIA, LFA, IHC, WB   |
| <b>Clonality</b>         | Polyclonal  |
| <b>Expiration Date</b>   | 12 months from date of receipt.   |