



Product Datasheet

GIT1 Antibody (orb99082)

Description	Mouse Anti-Rat GIT1 Monoclonal IgG1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Immunogen	Fusion protein amino acids 375-770 (C-terminus) of rat GIT1
Target	GIT1
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	-20°C
Note	For research use only
Application notes	1 µg/ml of SMC-413 was sufficient for detection of GIT1 in 10 µg of rat brain lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG1
Clonality	Monoclonal
MW	90kDa
Uniprot ID	Q9Z272
NCBI	NP_114002.1
Entrez	83709
Dilution Range	WB (1:1000)
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of neuroblas...



Western Blot analysis of rat brain membr...