


Product Datasheet

HDAC3 (phospho-Ser424) antibody (orb99117)

| | |
|----------------------------|---|
| Description | HDAC3 (phospho-Ser424) antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ICC, IF, IHC-P, WB |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Polyclonal |
| Dilution Range | WB 1:500-1:2000, IHC-P 1:50-1:200, IF/ICC 1:100-1:500 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of NIH-3T3 whole c...