


Product Datasheet

Tyrosine Hydroxylase (phospho-Ser19) antibody (orb99127)

Description	Tyrosine Hydroxylase (phospho-Ser19) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat cell lys...