


Product Datasheet

IRF-3 (phospho-Ser385) antibody (orb99164)

Description	IRF-3 (phospho-Ser385) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB 1:500-1:2000, IHC-P 1:50-1:200, IF/ICC 1:100-1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HT29 cell lysat...