


Product Datasheet

Mouse IL-33 ELISA Kit (orb99231)

| | |
|----------------------------|---|
| Description | Mouse IL-33 ELISA Kit |
| Reactivity | Mouse |
| Range | 31.2 pg/ml - 2,000 pg/ml |
| Tested Applications | ELISA |
| Immunogen | Expression system for standard: E.coli; Immunogen sequence: S109-I266 |
| Storage | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| Note | For research use only |
| Application notes | Dilution ratio of 1:1, concentration in serum and plasma is 0-260 pg/ml |
| Sample Types | cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA). |
| MW | 29991 MW |
| Uniprot ID | Q8BVZ5 |
| Sensitivity | 10 pg/ml |
| Expiration Date | Please enquire. |



Mouse IL-33
PicoKine
ELISA Kit
standard ...