


Product Datasheet

JAK1(Phospho-Tyr1022) antibody (orb99251)

Description	Rabbit polyclonal to JAK1. Janus kinase 1 (JAK1), is a member of a new class of protein-tyrosine kina...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000,IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis on MCF7 cell lysat...