


Product Datasheet

FOXO1/3/4-pan (Phospho-Thr24/32) antibody (orb99301)

Description	Rabbit polyclonal to FOXO1-3-4-pan.FOXO4 transcription factor AFX1 containing 1 fork-head domain.May...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000,IHC-P:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis on 293 treated wit...