


Product Datasheet

HDAC5 (Phospho-Ser259) antibody (orb99310)

Description	Rabbit polyclonal to HDAC5.HDAC9 a transcriptional regulator of the histone deacetylase family, subf...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis on HepG2 using pho...