



Product Datasheet

BCL2 antibody (orb99416)


Description	Rabbit polyclonal antibody to BCL2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB 1:500-1:2000, IHC-P 1:50-1:200, IF/ICC 1:100-1:500, ELISA(peptide): 1:20000-1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis of EC304 cells(Lan...



Western blot analysis of extracts from M...



Western blot analysis of 293 cells treat...