



Product Datasheet

mTOR antibody (orb99435)


Description	Rabbit polyclonal antibody to mTOR
Species/Host	Rabbit
Reactivity	Fish, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF/ICC 1:100-1:500, ELISA: 1:20000-1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa whole cell...



Western blot analysis of Mouse spleen us...



Immunohistochemical staining of human pr...