



## Product Datasheet

### GATA4 antibody (orb99437)


<b>Description</b>	GATA4 antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, ICC, IF, WB
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Dilution Range</b>	WB 1:500-1:2000, IF/ICC 1:100-1:500, ELISA: 1:20000-1:40000
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of HeLa whole cell...



Immunofluorescence analysis of 293 cells...



Immunofluorescence analysis of Hep G2 ce...