



Product Datasheet

SOX9 antibody (orb99439)


Description	SOX9 antibody
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, IP, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB 1:500-1:2000, IHC-P 1:50-1:200, IF/ICC 1:100-1:500, IP: 1:200, ELISA: 1:20000-1:40000
Expiration Date	12 months from date of receipt.



Western blot analysis of COS-7 treated w...



Immunohistochemical analysis of human br...



Immunofluorescence analysis of 293 cells...